

Created: 20200210
Last Modified: 20200212

The images acquired from the llsm needs to be preprocessed. Currently, we use the code “preprocessing_perAcquisition_ver3.m”, which will do the preprocessing when all images have completed the acquisition (in a single Ex folder). To carry this out successfully, we will need to correctly set up the folder structure.

Folder structure for the LLSM:

- D:
 - {[Year][Month][Day]_[InnerNA]_[OuterNA]_[Camera]_[Biologist]}
 - LLSCalibrations
 - illum
 - sB_noDither
 - {[sBnD_...].tif}
 - nB_noDither
 - {[nBnD_...].tif}
 - nB_Dither
 - {[nBD_...].tif}
 - chroma
 - {Ex[%02d]chroma[_{wavelength}nm_{power}mW]_[s-plane_spacing]}
 - {[chroma_...].tif}
 - bk
 - [bk_..._CamA_...].tif
 - [bk_..._CamB_...].tif
 - PSF
 - {Ex[%02d]psf[_{wavelength}nm_{power}mW]_[s-plane_spacing]}
 - {psf_..._.tif}
 - XZPSF
 - {[wavelength]XZPSF}
 - {CS[%1d]_[sampleName]_[sampleInfo]}
 - {Ex[%02d][keyword][_{wavelength}nm_{power}mW]_[s-plane_spacing]}
 - {[ex[%02d]_..._.tif}

Example

- D:
 - 20210210_p5_p55_sCMOS_Gustavo
 - LLSCalibrations
 - illum
 - sB_noDither
 - sBnD_..._CamA_..._488nm_... .tif
 - sBnD_..._CamB_..._560nm_... .tif
 - sBnD_..._CamA_..._642nm_... .tif
 - nB_noDither
 - sBnD_..._CamA_..._488nm_... .tif
 - sBnD_..._CamB_..._560nm_... .tif
 - sBnD_..._CamA_..._642nm_... .tif
 - nB_Dither
 - sBnD_..._CamA_..._488nm_... .tif
 - sBnD_..._CamB_..._560nm_... .tif
 - sBnD_..._CamA_..._642nm_... .tif
 - chroma
 - Ex00chroma_488nm_300mW_560nm_500mW_642nm_500mW_z0p2
 - chroma_...._488nm_... .tif
 - chroma_...._560nm_... .tif
 - chroma_...._642nm_... .tif
 - bk
 - bk_..._CamA_... .tif
 - bk_..._CamB_... .tif
 - PSF
 - Ex488psf1_488nm_300mW_z0p2
 - psf_...488nm_... .tif
 - Ex560psf1_488nm_300mW_z0p2
 - psf_...560nm_... .tif
 - Ex642psf1_488nm_300mW_z0p2
 - psf_...642nm_... .tif
 - XZPSF
 - 488XZPSF
 - 560XZPSF
 - 642XZPSF
 - CS1_SVGA_eGFP_NUP133_JF649_NUP107
 - Ex01_488nm_300mW_560nm_500mW_z0p2
 - Ex02_488nm_300mW_560nm_500mW_z0p2
 - Ex03_488nm_300mW_560nm_500mW_z0p2

- CS2_SVGA_eGFP_NUP133_JF649_NUP205
 - Ex04_488nm_300mW_560nm_500mW_z0p2
 - Ex05_488nm_300mW_560nm_500mW_z0p2
 - Ex06_488nm_300mW_560nm_500mW_z0p2
-

How to read:

- {values}
 - Means that there are multiple DEFINED value within values.
 - For example, an experiment with the channels: 488, 560 and 642; powers: 300, 500, 500
 - {wavelength} = 488560642
 - {_wavelength}nm = _488nm_560nm_642nm
 - {_wavelength}nm_{power}mW = _488nm_300mW_560nm_500mW_642nm_500mW
- [variable]
 - Means it is a NON-DEFINED variable. We also use regular expressions
 - Source folder name
 - For example, if todays date is 2021, Feb 10th
 - [Year] = 2021
 - [Month] = 02
 - [Day] = 10
 - For example, if the inner NA is 0.5 and the outer is 0.55
 - [innerNA] = p5
 - [outerNA]= p55
 - Camera will be either sCMOS or EMCCD. Lets say we are using sCMOS
 - [Camera] = sCMOS
 - Biologist will be the name of the biologist, lets say its Gustavo, then
 - [Biologist] = Gustavo
 - Experiments
 - [%02d] means any two digit number, so can be 00, 10, 32, 99, ... etc
 - Ex[%02d] means it can any of, but not limited to: Ex00, Ex10, Ex32, Ex99,... etc
 - Of course, if you have more than 100 unique experiments, you can use Ex[%03d]_ format
 - Ex100_

- Ex999_
 - [keyword]
 - Your personal keyword. You don't have to use it. You can use it to distinguish experiments. For example
 - good
 - nup205
 - lowSignal
 - incubated30min
 - " (empty string)
 - Ex[%02d][keyword]_. THIS MUST BE UNIQUE
 - Ex01good_
 - Ex10lowSignal_
 - Ex15SNR10_
 - Ex43incubated30min_
 - Ex90_
 - Combination of Both
 - Source folder
 - Examples of
 {[Year][Month][Day]_[InnerNA]_[OuterNA]_[Camera]_[Biologist]}:
 - 20200210_p5_p55_sCMOS_Gustavo
 - 20201222_p6_p5_EMCCD_Alex
 - 20210123_p5_p55_sCMOS_Anwesha
-

Notes

- If any of [Year][Month][Day]_[InnerNA]_[OuterNA]_[Camera]_[Biologist] are different, create a new folder
 - Do not combine same day experiments into a single folder if you have different biologist
 - You can copy and paste the calibrations folder from one to another if you are using exactly the same condition (same filter, imaging media, etc..)
- Each Ex[%02d][keyword] should have a unique cell/sample to be imaged
 - For example, if you image the same sample 3 times in a different folder, then you can name your folder something like (this is an example):
 - Ex01_...._
 - Ex01second_ ...
 - Ex01third_ ...
 - For example, when we ran NUP experiments, we first took a volume scan, then bleached the cell, then ran a time-series acquisition. The folder names can be:
 - Ex01initialVol_

- Ex01bleach_
 - Ex01timeSeries_
- Each Ex[%02d][keyword]_ MUST BE UNIQUE
- If you acquire a volume scan (sample scan or objective scan) you must create a folder with “Ex ... “ and save the tif images there